



Archiving Solutions

Secure Slide and Cassette Storage





Importance of Slide and Cassette Archiving

Archiving is an essential process in the laboratory to ensure sample security throughout the entire histology workflow.

Epredia's high quality cabinets are ideal for the safe storage of microscope slides and histology cassettes.

- Safely store slides and cassettes with Colorstack cabinets
- Securely manage samples in a dust-free environment
- Variety of sizes to meet all laboratory needs




Slide Storage

Slide Cabinets and Drawers



- ### Slide Filing Cabinet
- 14 drawer slide filing cabinet
 - Store 6000 76 x 26mm microscope slides with these 14-Drawer Cabinets


Item Number	Description	EMEA	AMER	APAC
E43/12GREEN	14 Drawer Slide Filing Cabinet (Green)	✓	✓	✓
E43/12GREY	14 Drawer Slide Filing Cabinet (Grey)	✓	✓	✓
E43/12BLUE	14 Drawer Slide Filing Cabinet (Blue)	✓	✓	✓
E43/13GREY	7 Drawer Slide Filing Cabinet (Grey)	✓		
E43/13BLUE	7 Drawer Slide Filing Cabinet (Blue)	✓		
E4313GRAY	7 Drawer Slide Filing Cabinet (Grey)		✓	



- ### Slide File
- Holds 3 x 1" (75 x 26 mm) slides
 - Store up to 100 slides (one per slot), 200 slides (two per slot), or 400 slides (slide tray removed)
 - Transparent hinged cover
 - Numbered slots for easy identification
 - Stackable for space-efficient storage


Item Number	Description	EMEA	AMER	APAC
B1000700BL	100 Slide File Box (Blue)		✓	

- Key**
- ✓ Denotes product in geography
 - AMER** North America and South America
 - EMEA** Europe, Middle East, Africa
 - APAC** Asia, Australia



- ### Colorstack Cabinet
- Colorstack cabinets have been designed to color-code the storage of 76 x 26 mm slides
 - Cabinet and drawer bodies are supplied in a cream epoxy-coated metal
 - Top covers available to match
 - Fit standard base stands

Item Number	Description	EMEA	AMER	APAC
E43/12C/BLUE	14 Drawer Colorstack Cabinet (Blue)	✓	✓	
E43/12C/GREY	14 Drawer Colorstack Cabinet (Grey)	✓	✓	
E4312C/TAN	14 Drawer Colorstack Cabinet (Tan)		✓	
E4312C/ORANGE	14 Drawer Colorstack Cabinet (Orange)		✓	
E43/12C/GREEN	14 Drawer Colorstack Cabinet (Green)	✓	✓	
E43/12C/PINK	14 Drawer Colorstack Cabinet (Pink)	✓		
E43/12C/TAN	14 Drawer Colorstack Cabinet (Tan)	✓		



- ### Lids and Covers for Cabinets and Bases
- Hold up to 5000 standard 3 × 1 in
 - (75 × 26mm) microscope slides in the Stackable Slide Storage Cabinet
 - E43/15 style bases are compatible with E43/12 cabinets


Item Number	Description	EMEA	AMER	APAC
E43/14GREY	Top Cover for Cabinets (Grey)	✓		
E43/14BLUE	Top Cover for Cabinets (Blue)	✓		
E4314GREEN	Top Cover for Cabinets (Green)		✓	✓
E43/14GREEN	Top Cover for Cabinets (Green)	✓		✓
E43/15	Base Stand for Cabinets (Black)	✓	✓	✓
E4315/HDWB	Heavy Duty Wheel Base for Standard Cabinet		✓	
E43/15/HDWB	Heavy Duty Wheel Base for Standard Cabinet	✓		✓



- ### Enamel Slide Storage Cabinet
- Hold up to 5000 standard 3 × 1 in
 - (75 × 26mm) microscope slides in the Stackable Slide Storage Cabinet
 - Finished in beige baked enamel and complete with drawers
 - Cabinets lock together when stacked
 - Metal base (ordered separately) is mounted on four rubber feet to protect bench surfaces and floors

Item Number	Description	EMEA	AMER	APAC
141001	Base for Stackable Slide Storage Cabinet (Metal)	✓	✓	
141000	Stackable Slide Storage Cabinet with Drawers (Enamel)	✓	✓	

Slide Tracking and Storage



Spring Spacer

- Separate slides while drying with these spacers, available in several varieties. Slides fit into slide filing cabinet drawers
- If slides are stored before they are completely dry, the spring spacer will separate them to stop them sticking


Item Number	Description	EMEA	AMER	APAC
E43/16	Spring Spacer for E43/12 Cabinets	✓	✓	✓
E4321	Spring Space for E43/13 Cabinets		✓	



Slide Filing Location Slips

- These color strips are used to identify the location of a slide that has been removed from a cabinet

Item Number	Description	EMEA	AMER	APAC
E43/25	Slide Filing Location Slips (White)	✓	✓	
E43/26	Slide Filing Location Slips (Pink)	✓		



Foam Dividers


- Use these foam dividers with Epredia Slide Cabinets, Modular Filing System and Stackable Slide Storage Cabinets

Item Number	Description	EMEA	AMER	APAC
B1000152	Sponge Separator for 141000 Stackable Enamel Slide Drawers	✓	✓	

Horizontal Slide Storage



- Extremely durable plastic slide boxes to store 76 x 26 mm slides
- Three sizes available – 25, 50 or 100 slide capacity
- Lids separate on the 25 and 50 slide capacity boxes
- 50 capacity slide box has a plastic catch
- 100 capacity slide box has a metal clasp
- Three colors available – white, blue and red
- Cork lined with index card in lid



- 100 slide capacity
- Plastic construction
- Indexed to aid in retrieval of stored slides

Item Number	Description	EMEA	AMER	APAC
B1780	25 Position Slide Box (White)	✓	✓	✓
B1781	25 Position Slide Box (Red)	✓	✓	✓
B1782	25 Position Slide Box (Blue)	✓	✓	✓
B1783	50 Position Slide Box (White)	✓	✓	✓
B1784	50 Position Slide Box (Red)	✓	✓	
B1785	50 Position Slide Box (Blue)	✓	✓	
B1787	100 Position Slide Box (Red)	✓	✓	✓
B1788	100 Position Slide Box (Blue)	✓	✓	✓
58299	100 Position Slide Box (White)	✓	✓	✓

Slide Folders and Mailers



Plastic Slide Folder

- Plastic slide folders available for the easy color-coding of 76 x 26 mm slides
- 20 slide capacity
- Available in five colors
- Stackable with transparent cover
- 10 folders per case

Item Number	Description	EMEA	AMER	APAC
B1001200BL	20 Position Plastic Slide Folder (Blue)	✓	✓	✓
B1001200GR	20 Position Plastic Slide Folder (Green)	✓	✓	✓
B1001200WH	20 Position Plastic Slide Folder (White)	✓	✓	✓
B1001200YW	20 Position Plastic Slide Folder (Yellow)	✓	✓	✓
B1001200GY	20 Position Plastic Slide Folder (Grey)		✓	✓



Plastic Slide Case

- Plastic mailer holds five slides in a single-piece design with built in slide separators and a hinged, tight-fitting lid

Item Number	Description	EMEA	AMER	APAC
125	5 Position Plastic Slide Case	✓		



Slide Folders and Trays

- Slide Folders are made of heavy, reinforced cardboard

Item Number	Description	EMEA	AMER	APAC
E32	2 Position Cardboard Slide Folder	✓		
E49	9 Position Cardboard Slide Tray	✓		
E6/3	Slide Tray Box	✓		
1604	20 Position Cardboard Slide Folder	✓	✓	✓
2500	5 Position Cardboard Slide Folder		✓	

Cassette Storage



Blockstore Cabinet

- 3000 cassettes capacity (48 x 48 x 27.5cm) metal
- Permanently store up to 165 embedding rings or 250 cassettes per drawer with modular cassette storage drawers

Item Number	Description	EMEA	AMER	APAC
E43/22GREY	Blockstore Cassette Cabinet (Grey)	✓	✓	✓
E43/22BLUE	Blockstore Cassette Cabinet (Blue)	✓	✓	✓
E43/22GREEN	Blockstore Cassette Cabinet (Green)	✓		



Plastic Block Filing Cabinets

- Simple filing of up to 1500 blocks
- Plastic Block Filing Cabinets feature smooth sliding action and rigid construction

Item Number	Description	EMEA	AMER	APAC
B1001304	Cassette Storage Drawer	✓	✓	

Microscope Slides



Adhesion Slides

Ref	Description
J1800AMNZ	Superfrost™ Plus adhesion slides, White, 1440 per case
J7800AMNZ	Superfrost™ Plus adhesion slides, White with four clipped corners, 1440 per case
J1810AMNZ	Superfrost™ Plus adhesion slides, Blue, 1440 per case
J1820AMNZ	Superfrost™ Plus adhesion slides, Pink, 1440 per case
J1830AMNZ	Superfrost™ Plus adhesion slides, Yellow, 1440 per case
J1840AMNZ	Superfrost™ Plus adhesion slides, Green, 1440 per case
J2800AMNZ	Polysine™ adhesion slides, White, 1440 per case

Non-Adhesion Slides

Ref	Description
AG00008032E01MNZ50	Superfrost, White, 1000 per case
AF00008032E01MNZ50	Superfrost, White with four clipped corners, 1000 per case
AG00008132E01MNZ50	Superfrost, Blue, 1000 per case
AG00008232E01MNZ50	Superfrost, Pink, 1000 per case
AG00008332E01MNZ50	Superfrost, Yellow, 1000 per case
AG00008432E01MNZ50	Superfrost, Green, 1000 per case
AG00008532E01MNZ50	Superfrost, Orange, 1000 per case
AG00008632E01MNZ50	Superfrost, Violet, 1000 per case

SlideMate Slides

Ref	Description
TT-40418217-W	SlideMate White
TT-40418218-PS-W	SlideMate Plus White Adhesion




Cassettes



Laser Optimised Cassettes

Ref	Description
110-xx-10LM	Slotted cassette, loose lids, 1000 per case
110-xx-10LMI	Slotted cassette, loose lids, 1000 per case
115-xx-10LM	Slotted cassette, loose lids, 1000 per case
115-xx-10LMI	Slotted cassette, loose lids, 1000 per case
215-xx-10LM	Small square cassette, loose lids, 1000 per case
215-xx-10LMI	Small square cassette, loose lids, 1000 per case
315-xx-18LM	Round holes cassette, compatible with metal lid (FTML315), 1800 per case
315-xx-18LMI	Round holes cassette, compatible with metal lid (FTML315), 1800 per case
615-xx-10LMI	Small holes cassette, loose lids, 1000 per case
815-xx-18LM	Slotted cassette, attached lids, 1800 per case
815-xx-18LMI	Slotted cassette, attached lids, 1800 per case

Simply select your cassette and configuration number above, and add the two digit color code.



Example: Orange Embedding Cassette 315-06-18LM

Traditional Cassettes

Ref	Description
B851010xx	PrintMate™ embedding cassette, 10x75 tubes, no lids
B851040xx	PrintMate embedding cassette, 1000 loose, no lids
B851110xx	PrintMate biopsy cassette, 10x75 tubes, loose lids
B851120xx	PrintMate slotted cassette, 10x75 tubes, loose lids
B851729xx	Slotted cassette 2, 10x75 tubes, loose lids
B851734xx	Fine pore cassette, 10x75 tubes, attached lids
B851736xx	Thinline slotted cassette, 10x75 tubes
1001003	Metal lids for cassettes

Select your cassette and configuration number below, add the two digit colour code.
Example: Green EpreDia Slotted Cassette in Tubes = B851729GN



*full list of cassettes available on [epredia.com](https://www.epredia.com)





Contact your Epredia team

UK

sales.uk@epredia.com
service.uk@epredia.com

+44 (0)1928 628350
+44 (0)203 048 3806

Germany

sales.de@epredia.com
service.de@epredia.com

+49 (0) 214 5068 0026
+49 (0) 610 3599 6795

Italy

sales.italia@epredia.com
service.italia@epredia.com

+39 020 0702168
+39 020 0702168

Netherlands

sales.bnl@epredia.com
service.bnl@epredia.com

+31 (0) 76 799 6068
+31 (0) 76 799 6028

Belgium/Luxembourg

sales.bnl@epredia.com
service.bnl@epredia.com

+32 (0) 800 260 57
+32 (0) 800 260 57

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2023 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. M410089_EMEA_A4 R052023

epredia Enhancing precision
cancer diagnostics

POWERED BY SHANDON MICRON MENZEL-GLÄSER