



# Enhance traceability

Comprehensive sample identification demands advanced printing capabilities.

## **Epredia™ VEGA laser cassette printers offer:**

- High-resolution, durable printing
- Up to 100-character 2D barcodes
- Fast, six seconds per cassette printing
- Intelligent cassette selection from six-hopper carousel
- Quiet printing, < 65 dB





- High-resolution, 600-dpi barcode to include more information
- High-quality barcode printing may reduce instances of reprints and non-scannable barcodes



- LIS connectivity to streamline workflow
- Easy-to-use connection for smooth installation



- Reduce operational cost with no printhead replacements
- Robust laser technology consistently prints within six seconds for a faster turnaround time

## Ordering Information

	Epredia Ref.
<b>Epredia VEGA Laser Cassette Printer</b>	
Epredia VEGA laser cassette printer with all-in-one PC and six hoppers	VEGA05-FT-PCDI
Epredia VEGA laser cassette printer with 6 hoppers	8000038

<b>Accessories</b>	
Side delivery system for VEGA	9001140
Cassette sorter V2 for VEGA	8000042
Cassette collector for VEGA sorter, 4 x 25 cassette capacity	9001387
Duo unit V2 for VEGA	8000043
Trays couple for VEGA Duo, 2 x 10 cassette capacity	9001547
Hopper carousel with 6 hoppers, 480 cassette capacity	9001327
13.3 in. All-in-one PC	B8131005MD



For more information, contact your local Epredia Account Executive or visit [epredia.com](https://www.epredia.com)

