

# Efficient, Versatile, Ergonomic

**Epredia Histostar Embedding Center** 



# Epredia HistoStar Embedding Workstation

Cool, comfortable and in control



- Built-in literature tray for report storage
- Modular cold plate with 72-block capacity allows positioning to either the left or right of the embedding center for optimal workflow
- Large 5 L paraffin tank features a filtration system designed to prevent debris from entering the dispensing nozzle
- Convenient drawer for easily accessible storage

#### Engineered for efficiency

- Large 5 L paraffin capacity, a cold plate sized to fit 72 blocks capacity cold plate and storage for 300 cassettes.
- The Epredia<sup>™</sup> HistoStar<sup>™</sup>
   embedding workstation is
   designed to minimize service
   requirements, ensuring
   dependable performance for
   years to come.

#### Flexibility to fit your workflow

- The HistoStar workstation's adjustable paraffin-dispense paddle allows users to easily adjust the paddle to the most ergonomic position.
- A variety of nozzle lengths to fit for your workflow. Even Supa Mega cassettes can be embedded with ease.
- Integrated Epredia™ ParaTrimmer™ offers space-saving convenience.
   No more wax scraping at the microtome.
- The system's efficient wax residue pathway helps maintain a clean workspace.

#### Cool comfort

- Smooth, contoured form factor to eliminate uncomfortable pressure points.
- All user contact points are insulated to provide all-day comfort with a cool workspace.
- LED lighting uniformly illuminates
  the workspace reduces fatigue
  and minimizes errors without the
  clutter of awkward remote lamps.
  Lighting intensity is easily adjustable with five levels of illumination
  that can be individually controlled
  for the specimen and accessory
  areas.













- User-adjustable cold plate temperature designed to reduce the occurrence of cracked blocks
- High-capacity storage for 200 standard base molds,
   300 cassettes, 6 Epredia Excelsior™ or 4 Revos™ baskets
- LED lighting with five brightness settings offers optimal illumination for both the specimen and accessory area
- Accessory area offers ten tool positions with five different hole sizes
- ParaTrimmer design, optimized for maximum productivity
- Intuitive, multi-lingual, touch screen control system

#### **Epredia HistoStar Embedding Workstation Specifications** Embedding module (W, D, H) 25.6 x 23.6 x 16.1 in (650 x 600 x 410 mm) Dimensions Cold plate (W, D, H) 17 x 23.6 x 16.1 in (430 x 600 x 410 mm) Weight Embedding module 55.12 lbs (25 kg) Cold plate 44.10 lbs (20 kg) Operating noise level <65 dB Capacity, wax reservoir 1.32 gallon (5 litres) 122 - 158 °F (50 - 70 °C) Temperatures Wax reservoir Cold spot 41 °F (5 °C) Hot spot 122 - 158 °F (50 - 70 °C) Tissue storage 122 - 158 °F (50 - 70 °C) Mold storage 122 - 158 °F (50 - 70 °C) Cold plate Adjustable, 10 - 27 °F (-12 to -3 °C) with a tolerance of ±1.8 °F (1 °C)



#### **Ordering Information**

Description		Order Number
HistoStar dispensing	220 - 240 Vac 50/60Hz	A81000001
unit with cold plate	100 - 120 Vac 50/60Hz	A81000002
	220 - 240 Vac 50/60Hz (China Only)	A81000003
HistoStar cold plate	220 - 240 Vac 50/60Hz	A81000004
	100 - 120 Vac 50/60Hz	A81000005
	220 - 240Vac 50/60Hz (China Only)	A81000006
HistoStar embedding module	110 - 120 Vac / 220 - 240 Vac 50/60Hz	A81000007
	220 - 240 Vac 50/60Hz (China Only)	A81000008
Tamping tools kit, set of 2: 1 X 20 mm and 1 X 12 mm		A81010030
Spatula kit, set of two		A81010031
Cassette lid storage cover		A81010026
Magnifier assembly		A81010032
Footswitch		AP16629
Electric, heated forceps, red, 1 mm tip		A81010033
Electric, heated forceps, yellow, 2 mm tip		A81010034
Electric, heated forceps, blue, 4 mm tip		A81010035
Document holder kit		A81010036

### Efficient embedding requires an outstanding paraffin

Delivery, reliability, unmatched technical support and superior product quality are all part of what you receive with every case of Epredia™ Histoplast™ paraffin. Made of paraffin and plastic polymers, these three paraffins deliver excellent processing and yield thin, flat, wrinkle-free sections.



Contact your local Epredia representative for a side-by-side comparison.

Description	Quantity	Order Number
Histoplast PE	8 x 1 kg bags	8330
Histoplast IM	8 x 1 kg bags	8331
Histoplast LP	8 x 1 kg bags	8332

## Find out more at **epredia.com**

