



## CTM6 Coverslipper

### Efficiency for Time Savings

- Coverslip hundreds of slides per hour with a self-cleaning dispensing needle
- Quick and easy change of different sizes and styles of coverslips for versatility with different slides and tissue types

### User Safety

- load and unload rack without opening the hood
- Internal vapor extraction system with fan and permanent air-flow avoids any leakage of vapors at openings
- Integrated charcoal filtration
- Tool-free cleaning and maintenance

### Compact Form Factor

- Designed to fit almost any lab space and under a 30 in. or larger hood

Description	Order Number
CTM6 coverslipper	970010

Accessories	Order Number
Shallow tray for loading baskets	A79210092
Drying stand for five baskets	A79210056
Filter	9990610
Extraction ventilation kit	A79210080
Soft hair brush	P112257
Suction cups (3 / pk)	A79210081
Charcoal filter	184710

### Physical Specifications

Dimensions (D x W x H):  
17.7 x 18.1 x 14 in.  
(45 x 46 x 35.5 cm)

Weight:  
44.1 lbs / 20 kg

Power requirements:  
100-240 V; 50/60 Hz

### You'll need:

