



Meet the new face of performance

Tissue Cassettes from Epredia



The entire range of cassettes from EpreDia™ has been redesigned for better performance – from printing through laboratory processes to archiving

Printing

- The all-new surface finish and reinforced print faces are engineered to optimize print quality when used with PrintMate™ AS cassette printers from EpreDia
- All-new PrintMate AS-compatible tubes are precisely engineered to reduce jamming and prevent downtime

Grossing, Embedding and Processing

- Opening and closing forces have been carefully adjusted for ease of use, while safeguarding valuable specimens during processing
- Design and flow characteristics tested to ensure efficient processing

Sectioning

- The new cassettes are precision engineered for a solid connection to EpreDia microtomes

Archiving


- The new surface texture is designed to optimize scanability when used with the Arcos™ tissue block management system from EpreDia



The mark of performance

It's easy to identify cassettes designed for performance. Just look for the unique PrintMate Certified icon on the tab.





PrintMate Slotted Cassette from EpreDia

- Pore style: Slotted
- For routine specimens
- Lid hinge: Reversible

B851090 1,000 cassettes, loose, with loose-packed lids

B851120 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



Thinline Slotted Cassette

- Pore style: Slotted
- For routine specimens
- Lid hinge: Front edge

B851736 750 cassettes, 10 tubes of 75 cassettes, with color-matched, attached lids



PrintMate Biopsy Cassette from EpreDia

- Pore style: Small squares
- For biopsy specimens
- Lid hinge: Reversible

B851080 1,000 cassettes, loose, with loose-packed lids

B851110 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



Thinline Biopsy Cassette

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Front edge

B851737 750 cassettes, 10 tubes of 75 cassettes, with color-matched, attached lids



PrintMate Embedding Cassette from EpreDia

- Pore style: Round holes
- For routine specimens
- Lid hinge: Reversible

B851040 1,000 cassettes, loose, with no lids – for use with metal lids

B851010 750 cassettes, 10 tubes of 75 cassettes, with no lids – for use with metal lids

B851070 1,000 cassettes, loose, with loose-packed lids

B851100 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids

1001003 Metal lids for cassettes





**Slotted
Cassette 2**

- Pore style: Slotted
- For routine specimens
- Lid hinge: Front edge

B851739 1,000 cassettes, loose, with loose-packed lids
B851729 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



**TissueLoc™ Cassette
(Cassette IV) from EpreDia**

- Pore style: Slotted
- For routine specimens
- Lid hinge: Back edge
- Large tab for easy opening

B851954 1,000 cassettes, loose, with loose-packed lids
B851720 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



**HistoScreen™
Cassette from EpreDia**

- Pore style: Mesh
- For very small biopsy specimens
- Eliminates need for sponges or bags
- Lid hinge: Back edge

B851000 1,000 cassettes, loose, with loose-packed lids
B851001 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



**Biopsy
Cassette 2**

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Front edge

B851738 1,000 cassettes, loose, with loose-packed lids
B851730 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



**TissueLoc™ Biopsy
Cassette (Cassette IV)
from EpreDia**

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Back edge
- Large tab for easy opening

B851925 1,000 cassettes, loose, with loose-packed lids
B851750 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



Fine Pore Cassette

- Pore style: Small holes
- For very small biopsy specimens
- Lid hinge: Front edge

B851735 1,000 cassettes, loose, with color-matched, attached lids
B851734 750 cassettes, 10 tubes of 75 cassettes, with color-matched, attached lids

Ordering is easy.

Simply select your cassette and configuration number above, add the two digit color code from below.

Example: Light Green TissueLoc Biopsy Cassette in Tubes = B851925LG



WH
White




GY
Grey



TN
Tan



BE
Beige



YW
Yellow



GN
Green



LG
Light Green



LA
Lavender



BL
Blue



LI
Lilac



PK
Pink



OR
Orange



PrintMate AS on-demand cassette printer

Our cassettes are precisely engineered for optimal performance when paired with the PrintMate™ AS cassette printer from EpreDia.

Eliminate waste and streamline your labeling process by printing cassettes on-demand with the PrintMate AS at the grossing station.

- Save time – prints each cassette in ten seconds or less
- Print only the cassettes you need, right at the point of use
- Print barcodes on cassettes to automate downstream processes
- Eliminate errors endemic to handwritten processes
- Suitable for both on-demand and batch printing
- Optional collection system for efficient and organized batch printing
- Free your workspace, the PrintMate printer's compact footprint fits easily and integrates into your laboratory's workflow



Richard-Allan Scientific

4481 Campus Drive
Kalamazoo, MI 49008
United States
+1 (800) 522-7270

Find out more at [epredia.com](https://www.epredia.com)

 Enhancing precision
cancer diagnostics

POWERED BY SHANDON MICROM MENZEL-GLÄSER 